

Quanti-Tray Sealer



Year of Purchase - 2013

Cost - Rs. 5.49 Lac

Easy:

- No media preparation.
- No pipetting.
- No dilutions (for counts up to 200 with Quanti-Tray; 2,419 with Quanti-Tray/2000; and up to 2,272 with Quanti-Tray/Legioalert).

Rapid

- Less than 1 minute of hands-on time.
- Results in 24 hours or less (7 days for Legiolert).
- No confirmation necessary.

Accurate

- Detects down to one organism per 100 mL.
- 95% confidence limits better than 5- or 10-tube [Most Probable Number \(MPN\)](#).
- 95% confidence limits better than or comparable to membrane filtration (MF).

Cost effective

- Minimal equipment to purchase and maintain.
- No glassware to purchase and clean.
- Minimal bench space needed.

Specifications

Features	Specifications
Dimensions	30cm H x 27cm T x 39cm Bl
Temperature Range	0–32°C
Electrical Power Supply	230V, 50Hz, 3A (Modell 99-10896-00)
Weight	16 kg